

## QCA922A

The instrument monitors both the resonant frequency and the resonant resistance which are also reflected on the two analog outputs.


It is possible to connect the QCM device to a Bio-Logic potentiostat using a DB9-8BNC cable (catalog n°: 092-22/1).



092-QCA-FC



SE-CL3



eQCM		Catalog n°
Quartz crystal microbalance 27 MHz kit		SE-QCA922A
<b>Content</b>		
EQCM 27 MHz main unit and cable		SE-QCA922A-00
Male BNC/BNC cable (2 pieces, length 1 m)		COR28100
<b>Options</b>		
Connector from QCA to potentiostat		092-22/1
Low flow peristaltic pump		EL-AV-008
<b>Cells</b>		
Holder is needed to get a full qcm or eqcm set-up		092-QCA-FC
		Dip cell
		SE-CL3
		Well cell (ptfe)
		SE-CL4
		Well cell (peek)
		SE-CL4PK
		Transparent well cell
		SE-CL5
		Flow cell (90 µl) (ptfe)
		SE-CL6
		Flow cell (90 µl) transparent (ptfe)
		SE-CL7
		Flow cell (90 µl) (peek)
		SE-CL6PK
<b>Resonators</b>		
5 Mhz		Gold electrode (25 pieces)
		SE-5AU
9 Mhz	Etching finish	Gold electrode (30 pieces)
		SE-9AU-E/1
		Gold electrode (100 pieces)
		SE-9AU-E/3
	Polish finish	Gold electrode (30 pieces)
		SE-9AU-P/1
		Gold electrode (100 pieces)
		SE-9AU-P/3
Standard finish	Resonator and lead wire	Aluminum electrode (25 pieces)
		SE-9AL
		Gold electrode (25 pieces)
		SE-9AU
		Graphite electrode (25 pieces)
		SE-9C
		Copper electrode (25 pieces)
		SE-9CU
		Molybdenum electrode (25 pieces)
		SE-9MO
Nickel electrode (25 pieces)		
SE-9NI		
		Platinum electrode (25 pieces)
		SE-9PT
		Stainless steel (sus304) electrode (25 pieces)
		SE-9SS
Separated lead wire		Gold electrode (25 pieces)
		SE-9AU-S
		Platinum electrode (25 pieces)
		SE-9PT-S
Mirror finish	Resonator and lead wire	Aluminum electrode (25 pieces)
		SE-9AL-M
		Gold electrode (25 pieces)
		SE-9AU-M
		Gold electrode (500 pieces)
		SE-9AU-M2
		Graphite electrode (25 pieces)
		SE-9C-M
		Copper electrode (25 pieces)
		SE-9CU-M
Ito electrode (25 pieces)		
SE-9IT-M		
Nickel electrode (25 pieces)		
SE-9NI-M		
Platinum electrode (25 pieces)		
SE-9PT-M		
Silicon electrode (25 pieces)		
SE-9SI-M		
Stainless steel (sus304) electrode (25 pieces)		
SE-9SS-M		
Titanium electrode (25 pieces)		
SE-9TI-M		
Separated lead wire		Gold electrode (25 pieces)
		SE-9AU-MS
		Ito electrode (25 pieces)
		SE-9IT-MS
		Platinum electrode (25 pieces)
		SE-9PT-MS
<b>Others</b>		
Resonators lead wire no sputter (50 pieces)		SE-LEAD/2
Resonators 9 mhz, no sputter	No lead wire (50 pieces)	SE-9Q/2
	Mirror finish no lead wire (50 pieces)	SE-9Q-M/2
Control software for qcm922, qcm922a (no potentiostat control)		SE-WQCM

Gold electrode



Platinum electrode




Stainless steel electrode



# Measurement examples

## Electropolymerization of pyrrol

The polypyrrol film was coated on an Au quartz (used as the working electrode) using cyclic voltammetry (20 cycles).



The quartz electrode was immersed in an acetonitrile solution ( $\text{Bu}_4\text{NPF}_6$  0.2 mol/L) containing a solution of 1 methylpyrrol monomer (0.01 mol/L).

## Polypyrrol film growth on the quartz working electrode


Fig. 1 presents polypyrrol film growth on the quartz electrode during successive cycles of cyclic voltammetry. The reversibility of the charge transfer in such a polymer film is often dependent on the deposition mode (quasi-reversible in this example). The growth is very regular but tends to slow down during the last cycles. This can be due to an interfacial depletion of the solution in methyl pyrrol species in the layer close to the electrode surface and to a saturation of the working electrode surface area.

- (1): G. Sauerbrey, Phys. Verh., 1957, 8, 113-114.
- (2): G. Sauerbrey, Z. Phys., 1959, 155, 206-222.
- (3): Application note #13. Section "Apps & literature of EC-Lab division".

## QCM measurements during the film growth


Fig. 1 shows the resonant frequency decrease and the resonant resistance increase while the polymer film is growing. Moreover the variation is dependent on the potential sweep resulting in a pseudo oscillation of frequency and resistance related to successive cycles. This plot can also be made versus potential (see fig. 1).

Fig. 1: overlaid frequency and current vs.  $E_{we}$  of the polymer film growth. Scanning at 100 mV/s between 0 and 1.018 V.




Here, it is interesting to note the frequency evolution versus potential on one cycle. The figure below shows the time evolution on one cycle.

Fig. 2: graphic zoom on one cycle showing the resonant frequency and the current density versus elapsed time<sup>(3)</sup>.



The mass calculation is done automatically in the EC-Lab® process data tool. More details can be found in the application note<sup>(3)</sup>.



## QCA922A specifications

Frequency range	5 MHz-30 MHz, resolution 0.01 Hz
Resonant resistance range	10-16 $\Omega$ , resolution 0.1 $\Omega$
$\Delta F$ output	Full scale: $\pm 10$ V (14 bit) $\pm 100$ Hz / $\pm 500$ kHz
Mass range	0.1 ng / Hz // 3.5 ng / Hz
Resistance output	Full scale: 0-10 V (14 bit) 10 $\Omega$ to 10 k $\Omega$
Gate time	Variable (10 ms / 20 ms / 100 ms / 1 s / 10 s)
Interface	USB

## Surface finishing

	Roughness	Electrode materials deposition
Standard finish	0.6 $\mu\text{m}$	Sputtered
Etching finish	0.6 $\mu\text{m}$	Vacuum deposition
Mirror finish	0.06 $\mu\text{m}$	Sputtered
Polish finish	0.06 $\mu\text{m}$	Vacuum deposition

To help you choose the type of finish, bear in mind that the lower the roughness, the larger the surface area.